

UTERINE SOUND

INSTRUCTIONS FOR USE

DESCRIPTION

The Uterine Sound is used to gauge the depth and position of the uterine cavity. It is also helpful with the placement of IUDs (Intra-uterine device) and the verification of the intrauterine position and length for HSG (Hysterosalpingography) procedures. The proximal portion of the sound is 3.5 mm in diameter tapering to 2.5 mm at the distal portion with a bulbous tip to reduce the risk of perforation.

The device is malleable to adapt to anatomical variations.

It is a sterile, disposable, single-use only device.

INDICATIONS

The Uterine Sound is used to gauge the depth and position of the uterine cavity. It is also helpful with the placement of IUDs and the verification of the intrauterine position and length for HSG procedures.

CONTRAINDICATIONS

The Uterine Sound should not be used when the following conditions are present:

- Infection, including but not limited to, Acute Cervicitis and Endometriosis.
- · Unexplained uterine bleeding.
- · A pregnancy or the suspicion of a pregnancy.

PRECAUTIONS

- · Check the packaging for damage which can compromise sterility.
- Visually inspect the sound for any surface damage prior to use.
- Disposable Uterine Sounds are designed for single-patient use only.
- This procedure should only be performed by a trained medical professional.
- The use of excessive force should be avoided.
- Minor bleeding, cervical inflammation, and abdominal cramps may occur.

WARNINGS

- Gently introduce and manipulate the tip of the sound through the cervical canal.
- · Forced entry could create a false tunnel and bleeding.

INSTRUCTIONS

The Uterine Sound is a single-use sound made of white, malleable, plastic material.

Prior to use, the sound can be curved between the thumb and finger in order to follow the curve of the cervical canal. The graduations are in centimeters and allow for the measurement of the uterine cavity.

Dispose of the Uterine Sound in accordance with all applicable Federal, State and local medical and hazardous waste practices.

www.medline.com Manufactured for: Medline Industries, Inc., Northfield, IL 60093 USA Made in Taiwan 1-800-MEDLINE RG17BTQ