

# QuickCheck Drug Screen Test Cup with Creatinine (Urine)

A rapid test for the simultaneous, qualitative detection of multiple drugs and drug metabolites in human urine.

## INTENDED USE

The QuickCheck Drug Screen Test Cup is a rapid chromatographic immunoassay for the qualitative detection of multiple drugs and drug metabolites in urine at the following cut-off concentrations:

Test	Analyte	Cut-off (ng/mL)
Amphetamine (AMP 1000)	d-Amphetamine	1,000
Barbiturates (BAR 300)	Secobarbital	300
Benzodiazepines (BZO 300)	Oxazepam	300
Buprenorphine (BUP 10)	Buprenorphine	10
Cocaine (COC 300)	Benzoylcegonine	300
Marijuana (THC 50)	11-nor-Δ9-THC-9 COOH	50
Methamphetamine (MET 1000)	d-Methamphetamine	1,000
Opiate (OPI 2000)	Morphine	2,000
Phencyclidine (PCP 25)	Phencyclidine	25
Oxycodone (OXY 100)	Oxycodone	100
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP 100)	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	100
Fentanyl (FYL 20)	Norfentanyl	20

Configurations of the QuickCheck Drug Screen Test Cup come with any combination of the above listed drug analytes with creatinine. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

## SUMMARY

The QuickCheck Drug Screen Test Cup is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes monoclonal antibodies to selectively detect elevated levels of specific drugs in urine.

### Amphetamine (AMP)

Amphetamine is a Schedule I controlled substance under Canada's Controlled Drugs and Substances Act available by prescription (Dexedrine®) and is also available on the illicit market. Amphetamines are a class of potent sympathomimetic agents with therapeutic applications. They are chemically related to the human body's natural catecholamines: epinephrine and norepinephrine. Acute higher doses lead to enhanced stimulation of the central nervous system (CNS) and induce euphoria, alertness, reduced appetite, and a sense of increased energy and power. Cardiovascular responses to amphetamines include increased blood pressure and cardiac arrhythmias. More acute responses produce anxiety, paranoia, hallucinations, and psychotic behavior. The effects of amphetamines generally last 2-4 hours following use and the drug has a half-life of 4-24 hours in the body. About 30% of amphetamines are excreted in the urine in unchanged form, with the remainder as hydroxylated and deaminated derivatives.

The QuickCheck Drug Screen Test Cup yields a positive result when the concentration of amphetamines in urine exceeds cutoff level.

### Barbiturates (BAR)

Barbiturates are CNS depressants. They are used therapeutically as sedatives, hypnotics, and anticonvulsants barbiturates are almost always taken orally as capsules or tablets. The effects resemble those of intoxication with alcohol. Chronic use of barbiturates leads to tolerance and physical dependence.

Short-acting barbiturates taken at 400 mg/day for 2-3 months can produce a clinically significant degree of physical dependence. Withdrawal symptoms experienced during periods of drug abstinence can be severe enough to cause death.

Only a small amount (less than 5%) of most barbiturates are excreted unaltered in the urine.

The approximate detection time limits for barbiturates are:

Short acting (e.g. Secobarbital)	100 mg PO (oral)	4.5 days
Long acting (e.g. Phenobarbital)	400 mg PO (oral)	7 days <sup>2</sup>

The QuickCheck Drug Screen Test Cup yields a positive result when the concentration of barbiturates in urine exceeds cutoff level.

### Benzodiazepines (BZO)

Benzodiazepines are medications that are frequently prescribed for the symptomatic treatment of anxiety and sleep disorders. They produce their effects via specific receptors involving a neurochemical called gamma aminobutyric acid (GABA). Because they are safer and more effective, benzodiazepines have replaced barbiturates in the treatment of both anxiety and insomnia. Benzodiazepines are also used as sedatives before some surgical and medical procedures, and for the treatment of seizure disorders and alcohol withdrawal.

Risk of physical dependence increases if benzodiazepines are taken regularly (e.g., daily) for more than a few months, especially at higher than normal doses. Stopping abruptly can bring on such symptoms as trouble sleeping, gastrointestinal upset, feeling

unwell, loss of appetite, sweating, trembling, weakness, anxiety and changes in perception.

Only trace amounts (less than 1%) of most benzodiazepines are excreted unaltered in the urine; most of the concentration in urine is conjugated drug. The detection period for benzodiazepines in urine is 3-7 days.

The QuickCheck Drug Screen Test Cup yields a positive result when the concentration of benzodiazepines in urine exceeds cutoff level.

### Buprenorphine (BUP)

Buprenorphine is a potent analgesic often used in the treatment of opioid addiction. The drug is sold under the trade names Subutex™, Buprenex™, Temgesic™ and Suboxone™, which contain Buprenorphine HCl alone or in combination with Naloxone HCl. Therapeutically, Buprenorphine is used as a substitution treatment for opioid addicts. Substitution treatment is a form of medical care offered to opiate addicts (primarily heroin addicts) based on a similar or identical substance to the drug normally used. In substitution therapy, Buprenorphine is as effective as Methadone but demonstrates a lower level of physical dependence. Concentrations of free Buprenorphine and Norbuprenorphine in urine may be less than 1 ng/ml after therapeutic administration, but can range up to 20 ng/ml in abuse situations. The plasma half-life of Buprenorphine is 2-4 hours.<sup>7</sup> While complete elimination of a single dose of the drug can take as long as 6 days, the window of detection for the parent drug in urine is thought to be approximately 3 days.

Substantial abuse of Buprenorphine has also been reported in many countries where various forms of the drug are available. The drug has been diverted from legitimate channels through theft, doctor shopping, and fraudulent prescriptions, and been abused via intravenous, sublingual, intranasal and inhalation routes.

The QuickCheck Drug Screen Test Cup yields a positive result when the Buprenorphine in urine exceeds cutoff level.

### Cocaine(COC)

Cocaine is a potent central nervous system stimulant and a local anesthetic. Initially, it brings about extreme energy and restlessness while gradually resulting in tremors, over-sensitivity and spasms. In large amounts, cocaine causes fever, unresponsiveness, difficulty in breathing and unconsciousness.

Cocaine is often self-administered by nasal inhalation, intravenous injection and free-base smoking. It is excreted in the urine in a short time primarily as benzoylecgonine.<sup>3,4</sup> Benzoylecgonine, a major metabolite of cocaine, has a longer biological half-life (5-8 hours) than cocaine (0.5-1.5 hours), and can generally be detected for 24-48 hours after cocaine exposure.<sup>4</sup>

The QuickCheck Drug Screen Test Cup yields a positive result when the concentration of benzoylecgonine in urine exceeds cutoff level.

### Marijuana (THC)

THC (Δ9-tetrahydrocannabinol) is the primary active ingredient in cannabis (marijuana). When smoked or orally administered, THC produces euphoric effects. Users have impaired short-term memory and slowed learning. They may also experience transient episodes of confusion and anxiety. Long-term, relatively heavy use may be associated with behavioral disorders. The peak effect of marijuana administered by smoking occurs in 20-30 minutes and the duration is 90-120 minutes after one cigarette. Elevated levels of urinary metabolites are found within hours of exposure and remain detectable for 3-10 days after smoking. The main metabolite excreted in the urine is 11-nor-Δ9-tetrahydrocannabinol-9-carboxylic acid (THC-COOH).

The QuickCheck Drug Screen Test Cup yields a positive result when the concentration of THC-COOH in urine exceeds cutoff level.

### Methamphetamine (MET)

Methamphetamine is an addictive stimulant drug that strongly activates certain systems in the brain. Methamphetamine is closely related chemically to Amphetamine, but the central nervous system effects of Methamphetamine are greater. Methamphetamine is made in illegal laboratories and has a high potential for abuse and dependence. The drug can be taken orally, injected, or inhaled. Acute higher doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, reduced appetite, and a sense of increased energy and power. Cardiovascular responses to Methamphetamine include increased blood pressure and cardiac arrhythmias. More acute responses produce anxiety, paranoia, hallucinations, psychotic behavior, and eventually, depression and exhaustion.

The effects of Methamphetamine generally last 2-4 hours and the drug has a half-life of 9-24 hours in the body. Methamphetamine is excreted in the urine primarily as Amphetamine, and oxidized and deaminated derivatives. However, 10-20% of Methamphetamine is excreted unchanged. Thus, the presence of the parent compound in the urine indicates Methamphetamine use. Methamphetamine is generally detectable in the urine for 3-5 days, depending on urine pH level.

The QuickCheck Drug Screen Test Cup yields a positive result when the methamphetamine in urine exceeds cutoff level.

### Morphine/Opiate (OPI)

The QuickCheck Drug Screen Test Cup yields a positive result when the concentration of morphine in urine exceeds 2,000 ng/mL. This is the suggested screening cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration (SAMHSA, USA).<sup>1</sup>

### Phencyclidine (PCP)

Phencyclidine, also known as PCP or Angel Dust, is a hallucinogen that was first marketed as a surgical anesthetic in the 1950's. It was removed from the market because patients receiving it became delirious and experienced hallucinations.

PCP is used in powder, capsule, and tablet form. The powder is either snorted or smoked after mixing it with marijuana or vegetable matter. PCP is most commonly administered by inhalation but can be used intravenously, intra-nasally, and orally. After low doses, the user thinks and acts swiftly and experiences mood swings from euphoria

to depression. Self-injurious behavior is one of the devastating effects of PCP.

PCP can be found in urine within 4 to 6 hours after use and will remain in urine for 7 to 14 days, depending on factors such as metabolic rate, user's age, weight, activity, and diet.<sup>6</sup> PCP is excreted in the urine as an unchanged drug (4% to 19%) and conjugated metabolites (25% to 30%).<sup>6</sup>

The QuickCheck Drug Screen Test Cup yields a positive result when the concentration of phencyclidine in urine exceeds 25 ng/mL. This is the suggested screening cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration (SAMHSA, USA).<sup>1</sup>

### Oxycodone (OXY)

Oxycodone is a semi-synthetic opioid with a structural similarity to codeine. The drug is manufactured by modifying thebaine, an alkaloid found in the opium poppy. Oxycodone, like all opiate agonists, provides pain relief by acting on opiate receptors in the spinal cord, brain, and possibly directly in the affected tissues. Oxycodone is prescribed for the relief of moderate to high pain under the well-known pharmaceutical trade names of OxyContin®, Tylox®, Percodan® and Percocet®. While Tylox®, Percodan® and Percocet® contain only small doses of oxycodone hydrochloride combined with other analgesics such as acetaminophen or aspirin, OxyContin consists solely of oxycodone hydrochloride in a time-release form. Oxycodone is known to metabolize by demethylation into oxymorphone and noroxycodone. In a 24-hour urine, 33-61% of a single, 5 mg oral dose is excreted with the primary constituents being unchanged drug (13-19%), conjugated drug (7-29%) and conjugated oxymorphone (13-14%). The window of detection for oxycodone in urine is expected to be similar to that of other opioids such as morphine.

The QuickCheck Drug Screen Test Cup yields a positive result when oxycodone in urine exceeds 100 ng/mL. This is the suggested screening cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration (SAMHSA, USA).<sup>1</sup>

### 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)

Methadone is an unusual drug in that its primary urinary metabolites (EDDP and EMDP) are cyclic in structure, making them very difficult to detect using immunoassays targeted to the native compound.<sup>10</sup> Exacerbating this problem, there is a subsection of the population classified as "extensive metabolizers" of methadone. In these individuals, a urine specimen may not contain enough parent methadone to yield a positive drug screen even if the individual is in compliance with their methadone maintenance. EDDP represents a better urine marker for methadone maintenance than unmetabolized methadone.

The QuickCheck Drug Screen Test Cup yields a positive result when the concentration of EDDP in urine exceeds cutoff level.

### Fentanyl (FYL)

Fentanyl is a potent, highly-addictive opioid narcotic that is used to treat severe chronic pain and is a μ special opiates receptor stimulant. Fentanyl is one of the varieties that have been listed in management of United Nations "Single Convention of narcotic drug in 1961". Among the opiates agents that under international control, fentanyl is one of the most commonly used to cure moderate to severe pain<sup>1</sup>. After continuous injection of fentanyl, the sufferer will have the performance of protracted opioid abstinence syndrome, such as ataxia and irritability etc.<sup>2,3</sup> which presents the addiction after taking fentanyl in a long time. Compared with drug addicts of amphetamine, drug addicts who take fentanyl mainly have got the possibility of higher infection rate of HIV, more dangerous injection behavior and more lifelong medication overdose <sup>4</sup>.

The QuickCheck Drug Screen Test Cup yields a positive result when the concentration of fentanyl in urine exceeds cutoff level.

## ADULTERATION

Adulteration is the tampering of a urine specimen with the intention of altering the test results. The use of adulterants can cause false negative results in drug tests by either interfering with the screening test and/or destroying the drugs present in the urine. Dilution may also be employed in an attempt to produce false negative drug test results. **Creatinine** is a waste product of creatine; an amino-acid contained in muscle tissue and found in urine.<sup>2</sup> A person may attempt to foil a test by drinking excessive amounts of water or diuretics such as herbal teas to "flush" the system. Creatinine is a way to check for dilution and flushing, which are the most common mechanisms used in an attempt to circumvent drug testing. Low creatinine levels may indicate dilute urine. Absence of creatinine (<5 mg/dl) is indicative of a specimen not consistent with human urine.

## PRINCIPLE

During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody. The antibody will then react with the drug-protein conjugate and a visible coloured line will show up in the test region of the specific drug dipstick. The presence of drug above the cut-off concentration will saturate all the binding sites of the antibody. Therefore, the coloured line will not form in the test region.

A drug-positive urine specimen will not generate a coloured line in the specific test region of the dipstick because of drug competition, while a drug-negative urine specimen will generate a line in the test region because of the absence of drug competition.

To serve as a procedural control, a coloured line will always appear at the control region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

## REAGENTS

Each test line contains anti-drug mouse monoclonal antibody and corresponding drug-protein conjugates. The control line contains goat anti-rabbit IgG polyclonal antibodies and rabbit IgG.

## CREATININE TEST REAGENT

Adulteration Pad	Reactive indicator	Buffers and non-reactive ingredients
Creatinine	0.04%	99.95%

## PRECAUTIONS

- For healthcare professionals including professionals at point of care sites.
- For in vitro diagnostic use only. The test cup should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test cup should be discarded according to federal, provincial and local regulations.

## STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2°C to 30°C. The test is stable through the expiration date printed on the sealed pouch. The test cups must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

## SPECIMEN COLLECTION AND PREPARATION

### Urine Assay

The user can urinate directly into the QuickCheck Test Cup. Alternatively, urine specimen can be collected in a separate clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing.

### Specimen Storage

Urine specimens in a separate container may be stored at 2°C to 8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed well before testing. When testing with creatinine strip, storage of urine specimens should not exceed 2 hours at room temperature or 4 hours refrigerated prior to testing.

Urine specimen in Test Cups should not be stored and must be discarded immediately after testing.

## MATERIALS PROVIDED

- Drug Screen Test Cups
- Instruction For Use

## MATERIALS REQUIRED BUT NOT PROVIDED

- Timer
- Security Seal

## DIRECTIONS FOR USE

Allow the test cup, urine specimen (if stored), and/or controls to reach room temperature (15°C to 30°C) prior to testing.

- Bring the pouch containing the test cup to room temperature before opening it. Remove the cup from the sealed pouch and use it within one hour.
- Label the cup with client name or identification (ID) and date.
- Donor provides urine specimen. The amount of urine must reach the "Minimum Urine Level" mark on the Test Cup.
- Technician replaces and secures cap while the cup is on a flat surface.
- For freshly voided urine, check the temperature indicator label on the Test Cup up to 4 minutes after specimen collection. A green color will appear to indicate the temperature of the urine specimen. The proper range for an unadulterated specimen is 32°C to 38°C.
- If needed, the technician dates and initials the security seal and attaches the security seal over the cup cap. (Security seal is available separately from Innovatek Medical Inc.)
- Technician peels off label to reveal adulteration strip(s), if applicable.
- Technician peels off the label on the multi-drug test cup to view results.
- Read the creatinine between 3-5 minutes with the help of colour chart below. Refer to your Drug Free Policy for guidelines on adulterated specimens. We recommend not to interpret the drug test results and either retest the urine or collect another specimen in case of any abnormal result for the creatinine strip.
- The drug strip result should be read at 5 minutes. Do not interpret the result after 10 minutes.

## INTERPRETATION OF RESULTS

Please refer to the illustration below

### NEGATIVE:\*

A coloured line appears in the Control region (C) and coloured line appear in the Test region (T). This negative result means that the concentration in the urine sample is below the designated cut-off level for a particular drug tested.

\*NOTE: The shade of the coloured lines in the Test region (T) may vary. The result should be considered negative whenever there is even a faint line.

### POSITIVE:


A coloured line appears in the Control region (C) and NO line appears in the Test region (T). The positive result means that the drug concentration in the urine sample is greater than the designated cut-off for a specific drug.

### INVALID:


No line appears in the Control region (C). Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for Control line failure. Read the directions again and repeat the test with a new test card. If the result is still invalid, contact your manufacturer.



Read the drug strips at 5 minutes



Interpret creatinine strip between 3-5 minutes



## INTERPRETATION OF CREATININE TEST

(Please refer to the colour chart above)

Semi Quantitative results are obtained by visually comparing the reacted colour blocks on the creatinine strip to the printed colour blocks on the colour chart. No instrumentation is required. If the test pad colour matches that under "Normal" it indicates the urine specimen has not been adulterated.

## QUALITY CONTROL

A procedural control is included in the test. A line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

## LIMITATIONS

- The QuickCheck Drug Screen Test Cup provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.<sup>1,10</sup>
- There is a possibility that technical or procedural errors, as well as interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result does not indicate level or intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- This test does not distinguish between drugs of abuse and certain medications.

## CREATININE TEST LIMITATIONS

Normal creatinine level is between 20 and 350 mg/dL. Under rare conditions, certain kidney diseases may show dilute urine.

## EXPECTED VALUES

The negative result indicates that the drug concentration is below the detectable level. Positive result means the concentration of drug is above the detectable level.

## PERFORMANCE CHARACTERISTICS

### Accuracy

A side-by-side comparison was conducted using the QuickCheck Drug Screen Test Cup and commercially available drug rapid tests. Testing was performed on approximately 250 specimens per drug type previously collected from subjects presenting for Drug Screen Testing. Presumptive positive results were confirmed by GC/MS.

Method	QuickCheck Drug Screen Test Cup	GC/MS		% agreement with GC/MS
		Positive	Negative	
AMP 1000	Positive	103	3	98.1%
	Negative	2	142	97.9%
BAR 300	Positive	98	2	96.1%
	Negative	4	146	98.6%
BZO 300	Positive	121	1	98.4%
	Negative	2	126	99.2%
BUP 10	Positive	105	0	99.1%
	Negative	1	144	>99.9%
COC 300	Positive	111	3	98.2%
	Negative	2	134	97.8%
THC	Positive	92	3	97.9%

Method	QuickCheck Drug Screen Test Cup	GC/MS		% agreement with GC/MS
		Positive	Negative	
50	Negative	2	153	98.1%
	Positive	94	3	97.9%
THC 30	Negative	2	151	98.1%
	Positive	76	5	96.2%
MET 1000	Negative	3	166	97.1%
	Positive	117	8	96.7%
OPI 2000	Negative	4	121	93.8%
	Positive	85	5	92.4%
PCP 25	Negative	7	153	96.8%
	Positive	84	1	97.7%
OXY 100	Negative	2	163	99.4%
	Positive	95	5	96.9%
EDDP 100	Negative	3	147	96.7%
	Positive	79	1	98.8%
FYL 20	Negative	1	169	99.4%
	Negative	2	64	97.0%

## % Agreement with Commercial Kit

	AMP 1000	BAR 300	BZO 300	BUP 10	COC 300	THC 50
Positive Agreement	>99.9% %	>99.9% %	>99.9% %	>99.9% %	>99.9% %	>99.9% %
Negative Agreement	>99.9% %	>99.9% %	>99.9% %	>99.9% %	>99.9% %	>99.9% %
Total Results	>99.9% %	>99.9% %	>99.9% %	>99.9% %	>99.9% %	>99.9% %

	MET 1000	OPI 2000	PCP 25	OXY 100	EDDP 100	FYL 20
Positive Agreement	>99.9% *	* %	>99.9% %	* %	* %	* %
Negative Agreement	>99.9% %	* %	>99.9% %	* %	* %	* %
Total Results	>99.9% %	* %	>99.9% %	* %	* %	* %

\* Note: Based on GC/MS data instead of Commercial Kit.

## Precision

A study was conducted at three hospitals by laypersons using three different lots of product to demonstrate the within run, between run and between operator precision. An identical card of coded specimens, containing drugs at concentrations of ± 50% and ± 25% cut-off level, was labeled, blinded and tested at each site. The results are given below:

### AMPHETAMINE (AMP 1000)

Amphetamine conc. (ng/mL)	n per site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	10	10	0	10	0	10	0
500	10	10	0	10	0	10	0
750	10	9	1	8	2	9	1
1,250	10	1	9	2	8	2	8
1,500	10	0	10	0	10	0	10

### BARBITURATES (BAR 300)

Secobarbital conc. (ng/mL)	n per site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	10	10	0	10	0	10	0
150	10	10	0	10	0	10	0
225	10	9	1	8	2	9	1
375	10	2	8	1	9	2	8
450	10	0	10	0	10	0	10

### BENZODIAZEPINES (BZO 300)

Oxazepam conc. (ng/mL)	n per site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	10	10	0	10	0	10	0
150	10	10	0	10	0	10	0
225	10	9	1	9	1	9	1
375	10	1	9	1	9	1	9
450	10	0	10	0	10	0	10

### BUPRENORPHINE (BUP 10)

Buprenorphine conc. (ng/mL)	n per site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	10	10	0	10	0	10	0
5	10	10	0	10	0	10	0
7.5	10	9	1	9	1	8	2
12.5	10	1	9	1	9	1	9
15	10	0	10	0	10	0	10

COCAINE (COC 300)								
Benzoylcocgonine conc. (ng/mL)	n per site	Site A		Site B		Site C		
		-	+	-	+	-	+	
0	10	10	0	10	0	10	0	0
150	10	10	0	10	0	10	0	0
225	10	9	1	9	1	9	1	9
375	10	1	9	1	9	1	9	9
450	10	0	10	0	10	0	10	0

MARIJUANA (THC 50)								
11-nor- $\Delta^9$ -COOH conc. (ng/mL)	n per site	Site A		Site B		Site C		
		-	+	-	+	-	+	
0	10	10	0	10	0	10	0	0
25	10	10	0	10	0	10	0	0
37.5	10	9	1	8	2	9	1	9
62.5	10	1	9	1	9	2	8	8
75	10	0	10	0	10	0	10	0

METHAMPHETAMINE (MET 1000)								
Methamphetamine conc. (ng/mL)	n per site	Site A		Site B		Site C		
		-	+	-	+	-	+	
0	10	10	0	10	0	10	0	0
500	10	10	0	10	0	10	0	0
750	10	9	1	9	1	9	1	9
1,250	10	1	9	2	8	1	9	9
1,500	10	0	10	0	10	0	10	0

MORPHINE/OPIATE (OPI 2000)								
Morphine conc. (ng/mL)	n per site	Site A		Site B		Site C		
		-	+	-	+	-	+	
0	10	10	0	10	0	10	0	0
1,000	10	10	0	10	0	10	0	0
1,500	10	9	1	9	1	9	1	9
2,500	10	1	9	1	9	1	9	9
3,000	10	0	10	0	10	0	10	0

PHENCYCLIDINE (PCP 25)								
Phencyclidine conc. (ng/mL)	n per site	Site A		Site B		Site C		
		-	+	-	+	-	+	
0	10	10	0	10	0	10	0	0
12.5	10	10	0	10	0	10	0	0
18.75	10	8	2	9	1	9	1	9
31.25	10	1	9	1	9	1	9	9
37.5	10	0	10	0	10	0	10	0

OXYCODONE (OXY 100)								
Oxycodone conc. (ng/mL)	n per site	Site A		Site B		Site C		
		-	+	-	+	-	+	
0	10	10	0	10	0	10	0	0
50	10	10	0	10	0	10	0	0
75	10	9	1	9	1	9	1	9
125	10	1	9	1	9	1	9	9
150	10	0	10	0	10	0	10	0

2-ETHYLIDENE-1,5-DIMETHYL-3,3-DIPHENYLPYRROLIDINE (EDDP 100)								
EDDP conc. (ng/mL)	n per site	Site A		Site B		Site C		
		-	+	-	+	-	+	
0	10	10	0	10	0	10	0	0
50	10	10	0	10	0	10	0	0
75	10	9	1	9	1	9	1	9
125	10	1	9	1	9	1	9	9
150	10	0	10	0	10	0	10	0

FENTANYL (FYL 20)								
FYL conc. (ng/mL)	n per site	Site A		Site B		Site C		
		-	+	-	+	-	+	
0	10	10	0	10	0	10	0	0
10	10	10	0	10	0	10	0	0
15	10	9	1	9	1	9	1	9
25	10	1	9	1	9	1	9	9
30	10	0	10	0	10	0	10	0

**Analytical Sensitivity**  
A drug-free urine pool was spiked with drugs at the listed concentrations. The results are summarized below.

Drug Concentration Cut-off Range	AMP 1000		BAR 300		BZO 300		BUP 10		COC 300		THC 50	
	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0
-50% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	26	4	27	3	27	3	26	4	26	4	26	4
Cut-off	15	15	16	14	15	15	14	16	13	17	14	16

+25% Cut-off	3	27	4	26	3	27	3	27	3	27	3	27
+50% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30
+300% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30

Drug Concentration Cut-off Range	MET 1000		OPI 2000		PCP 25		OXY 100		EDDP 100		FYL 20	
	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0
-50% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	27	3	27	3	25	5	27	3	26	4	27	3
Cut-off	16	14	14	16	15	15	15	15	15	15	14	16
+25% Cut-off	3	27	4	26	3	27	4	26	3	27	4	26
+50% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30
+300% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30

**Analytical Specificity**  
The following table lists the concentrations of compounds (ng/mL) that are detected as positive in urine by the QuickCheck Drug Screen Test Cup at 5 minutes.

Analytes	Concentration (ng/mL)	Analytes	Concentration (ng/mL)
<b>AMPHETAMINE (AMP 1000)</b>			
D,L-Amphetamine sulfate	300	Phentermine	1,000
L-Amphetamine	25,000	Maprotiline	50,000
(±) 3,4-Methylenedioxyamphetamine	500	Methoxyphenamine	6,000
		D-Amphetamine	1,000
<b>BARBITURATES (BAR 300)</b>			
Amobarbital	5,000	Alphenol	600
5,5-Diphenylhydantoin	8,000	Aprobarbital	500
Allobarbital	600	Butobarbital	200
Barbital	8,000	Butalbital	8,000
Talbutal	200	Butethal	500
Cyclopentobarbital	30,000	Phenobarbital	300
Pentobarbital	8,000	Secobarbital	300
<b>BENZODIAZEPINES (BZO 300)</b>			
Alprazolam	100	Bromazepam	900
a-hydroxyalprazolam	1,500	Chlordiazepoxide	900
Clobazam	200	Nitrazepam	200
Clonazepam	500	Norchlordiazepoxide	100
Clorazepatedipotassium	500	Nordiazepam	900
Delorazepam	900	Oxazepam	300
Desalkylflurazepam	200	Temazepam	100
Flunitrazepam	200	Diazepam	300
(±) Lorazepam	3,000	Estazolam	6,000
RS-Lorazepam-glucuronide	200	Triazolam	3,000
Midazolam	6,000		
<b>BUPRENORPHINE (BUP 10)</b>			
Buprenorphine	10	Norbuprenorphine	50
Buprenorphine 3-D-Glucuronide	50	Norbuprenorphine 3-D-Glucuronide	100
<b>COCAINE (COC 300)</b>			
Benzoylcocgonine	300	Cocacethylene	20,000
Cocaine HCl	200	Ecgonine	80,000
<b>MARIJUANA (THC 50)</b>			
Cannabinol	35,000	$\Delta^8$ -THC	17,000
11-nor- $\Delta^8$ -THC-9 COOH	30	$\Delta^9$ -THC	17,000
11-nor- $\Delta^9$ -THC-9 COOH	50		
<b>METHAMPHETAMINE (MET 1000)</b>			
p-Hydroxymethamphetamine	25,000	(±)-3,4-Methylenedioxy-methamphetamine	12,500
D-Methamphetamine	1,000		
L-Methamphetamine	20,000	Mephentermine	50,000
<b>MORPHINE/OPIATE (OPI 2000)</b>			
Codaine	2,000	Morphine	2,000
Ethylmorphine	3,000	Norcodeine	25,000
Hydrocodone	50,000	Normorphone	50,000
Hydromorphone	15,000	Oxycodone	25,000
Levorphanol	25,000	Oxymorphone	25,000
6-Monoacetylmorphine	3,000	Procaine	50,000
Morphine 3- $\beta$ -D-glucuronide	2,000	Thebaine	25,000
<b>PHENCYCLIDINE (PCP 25)</b>			
Phencyclidine	25	4-Hydroxyphencyclidine	12,500
<b>Oxycodone (OXY 100)</b>			
Oxycodone	100	Hydromorphone	50,000
Oxymorphone	300	Naloxone	25,000
Levorphanol	50,000	Naltrexone	25,000
Hydrocodone	25,000		

2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP 100)			
2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)			100
<b>Fentanyl (FYL 20)</b>			
Alfentanyl	600,000	Buspirone	15,000
Fenfuramine	50,000	Fentanyl	100
Norfentanyl	20	Sufentanyl	50,000

**Effect of Urinary Specific Gravity**  
Fifteen (15) urine samples of normal, high, and low specific gravity ranges (1.005-1.045) were spiked with drugs at 50% below and 50% above cut-off levels respectively. The QuickCheck Drug Screen Test Cup was tested in duplicate using fifteen drug-free urine and spiked urine samples. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

**Effect of Urinary pH**  
The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with drugs at 50% below and 50% above cut-off levels. The spiked, pH-adjusted urine was tested with the QuickCheck Drug Screen Test Cup. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

**Cross-Reactivity**  
A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or drug positive urine containing Amphetamine, Barbiturates, Benzodiazepines, Buprenorphine, Cocaine, Marijuana, Methamphetamine, Morphine, Phencyclidine, Oxycodone, EDDP, Fentanyl. The following compounds show no cross-reactivity when tested with the QuickCheck Drug Screen Test Cup at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds			
Acetophenetidin	Cortisone	Zomepirac	d-Pseudoephedrine
N-Acetylprocainamide	Creatinine	Ketoprofen	Quinidine
Acetylsalicylic acid	Deoxycorticosterone	Labeltol	Quinine
Aminopyrine	Dextromethorphan	Loperamide	Salicylic acid
Amoxicillin	Clofenac	Meprobamate	Serotonin
Ampicillin	Flunisal	Isoxsuprine	Sulfamethazine
l-Ascorbic acid	Digoxin	d,l-Propranolol	Sulindac
Apomorphine	Diphenhydramine	Nalidixic acid	Tetracycline
Aspartame	Ethyl-p-aminobenzoate	Naproxen	Tetrahydrocortisone-
Atropine	$\beta$ -Estradiol	Niacinamide	3-acetate
Benzoic acid	Estrone-3-sulfate	Nifedipine	Tetrahydrocortisone
Bilirubin	Erythromycin	Norethindrone	Tetrahydrozoline
Flunitrazepam	Fenoprofen	Noscapine	Thiamine
Flunitrazepam	Furosemide	d,l-Octopamine	Thioridazine
Caffeine	Genticic acid	Oxalic acid	d,l-Tyrosine
Cannabidiol	Hemoglobin	Oxolinic acid	Tolbutamide
Chloral hydrate	Hydralazine	Oxymetazoline	Triamterene
Chloramphenicol	Hydrochlorothiazide	Papaverine	Trifluoperazine
Chlorothiazide	Hydrocortisone	Penicillin-G	Trimethoprim
d,l-Chlorpheniramine	o-Hydroxyhippuric acid	Perphenazine	d,l-Tryptophan
Chlorpromazine	3-Hydroxytyramine	Phenelzine	Uric acid
Cholesterol	d,l-Isoproterenol	Prednisone	Verapamil
Clonidine			

**REFERENCES**

- Hawks RL, Chiang CN. *Urine Testing for Drugs of Abuse*. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.
- Tietz NW. *Textbook of Clinical Chemistry*. W.B. Saunders Company. 1986; 1735.
- Stewart DJ, Inaba T, Lucassen M, Kalow W. *Clin. Pharmacol. Ther.* April 1979; 25 ed: 464, 264-8.
- Ambre J. *J. Anal. Toxicol.* 1985; 9:241.
- Winger, Gail. *A Handbook of Drug and Alcohol Abuse*, Third Edition, Oxford Press, 1992, page 146.
- Robert DeCresce. *Drug Testing in the workplace*, 1989 page 114.
- Glass IB. *The International Handbook of Addiction Behavior*. Routledge Publishing, New York, NY. 1991; 216
- Cody JT. "Specimen Adulteration in drug urinalysis. *Forensic Sci. Rev.*, 1990, 2:63.
- Tsai SC. et al. *J. Anal. Toxicol.* 1998; 22 (6): 474
- Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, 6th Ed. Biomedical Publ., Foster City, CA 2002.
- Hardman JG, Limbird LE. *Goodman and Gilman's: The Pharmacological Basis for Therapeutics*. 10th Edition. McGraw Hill Medical Publishing, 2001; 208-209.
- Cumming E. (22 April 2010). "Mephedrone: Chemistry lessons". London: The Daily Telegraph. Retrieved 2010-09-14.
- "Drugs crackdown hailed a success". BBC News. 8 March 2010. Retrieved 2010-03-31.
- Kihara, Rhiannon; Day, Edward (May 2014). "Transient psychotic episodes following recreational use of NR-G-3". *Progress in Neurology and Psychiatry* 18 (3): 14–18. doi:10.1002/pnp.331. Retrieved 22 March 2015.
- Schifano F, Albanese A, Fergus S, Stair JL, Deluca P, Corazza O, Davey Z, Corkery J, Siemann H, Scherbaum N, Farre M, Torrens M, Demetrios Z, Ghodse AH. *Psychonaut Web, M.; Rednet Research, G.* (2010). "Mephedrone (4-methylmethcathinone; 'meow meow'): chemical, pharmacological and clinical issues". *Psychopharmacology* 214 (3):593–602. doi:10.1007/s00213-010-2070-

- x.ISSN 0033-3158. PMID 21072502.
16. "Assessment of Zopiclone". *World Health Organization. Essential Medicines and Health Products World Health Organization. p.9 (Section 5. Pharmacokinetics)*. Retrieved 5 December 2015.
  17. Kratzsch C, Tenberken O, Peters FT et al. Screening, library-assisted identification, and validated quantification of 23 benzodiazepines, flumazenil, zaleplone, zolpidem, and zopiclone in plasma by liquid chromatography/mass spectrometry with atmospheric pressure chemical ionization. *J. Mass Spec.* 39: 856-872, 2004.
  18. Gustavsen I, Al-Sammurraie M, Mørland J, Bramness JG. Impairment related to blood drug concentrations of zopiclone and zolpidem compared with alcohol in apprehended drivers. *Accid. Anal. Prev.* 41: 462-466, 2009.
  19. Baselt R, *Disposition of Toxic Drugs and Chemicals in Man*, 8th edition, Biomedical Publications, Foster City, CA, 2008, pp. 1677-1679.
  20. Calkins RF, Aktan GB, Hussain KL (1995). "Methcathinone: the next illicit stimulant epidemic?". *Journal of Psychoactive Drugs.* 27 (3): 277–85. doi:10.1080/02791072.1995.10472472. PMID 8594170.
  21. Methcathinone. <https://en.wikipedia.org/wiki/Methcathinone>.
  22. Bersani FS, Corazza O, Simonato P, Mylokosta A, Levani E, Lovaste R, Schifano F. (2013). "Drops of madness? Recreational misuse of tropicamide collyrium: early warning alerts from Russia and Italy". *General Hospital Psychiatry* 35 (5):571–3. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
  23. Malenka RC, Nestler EJ, Hyman SE (2009). "Chapter 15: Reinforcement and Addictive Disorders". In Sydor A, Brown RY. *Molecular Neuropharmacology: A Foundation for Clinical Neuroscience* (2nd ed.). New York: McGraw-Hill Medical. p. 375. ISBN 9780071481274.
  24. American Psychiatric Association (2013). "Substance-Related and Addictive Disorders". American Psychiatric Publishing. pp. 1–2. Retrieved 10 July 2015.
  25. Juliano LM, Griffiths RR (2004). "A critical review of caffeine withdrawal: empirical validation of symptoms and signs, incidence, severity, and associated features". *Psychopharmacology (Berl.)* 176 (1):1–29. doi:10.1007/s00213-004-2000-x. PMID 15448977. Archived from the original on 29 January 2012.
  26. Arnaud MJ. Pharmacokinetics and metabolism of natural methylxanthines in animal and man. *Handbook Exp Pharmacol* 2011; 200:33-91.
  27. Jeukendrup AE, Randell R. Fat burners: nutrition supplements that increase fat metabolism. *Obes Rev* 2011; 193:1-24.
  28. Cumming, E. (22 April 2010). "Mephedrone: Chemistry lessons". London: The Daily Telegraph. Retrieved 2010-09-14.
  29. "Drugs crackdown hailed a success". BBC News. 8 March 2010. Retrieved 2010-03-31.
  30. Kihara, Rhiannon; Day, Edward (May 2014). "Transient psychotic episodes following recreational use of NRG-3". *Progress in Neurology and Psychiatry* 18 (3): 14–18. doi:10.1002/pnp.331. Retrieved 22 March 2015.
  31. Schifano F, Albanese A, Fergus S, Stair JL, Deluca P, Corazza O, Davey Z, Corkery J, Siemann H, Scherbaum N, Farre M, Torrens M, Demetrovics Z, Ghodse AH, Psychonaut Web, M.; Rednet Research, G. (2010). "Mephedrone (4-methylmethcathinone; 'meow meow'): chemical, pharmacological and clinical issues". *Psychopharmacology* 214 (3):593–602. doi:10.1007/s00213-010-2070-x. ISSN 0033-3158. PMID 21072502.
  32. Work Group on Panic Disorder (January 2009). *APA Practice Guideline for the Treatment of Patients With Panic Disorder* (2nd ed.).
  33. "FDA approved labeling for Xanax revision 08/23/2011". Federal Drug Administration. 2011-08-23. p. 4. Retrieved 2011-09-14.
  34. "Xanax XR (Alprazolam) Clinical Pharmacology – Prescription Drugs and Medications". RxList. First Data Bank. July 2008.



## Innovatek Medical Inc.

3 - 1600 Derwent Way, Delta, BC, Canada. V3M 6M5



4215-DUA-1127A-125 QuickCheck 12 Test Cup  
 4215-DUA-1127A-150 QuickCheck 12 Test Cup  
 4215-DUA-1117A-125 QuickCheck 11 Test Cup  
 4215-DUA-1117A-150 QuickCheck 11 Test Cup  
 4215-DUA-197A-125 QuickCheck 9 Test Cup  
 4215-DUA-197A-150 QuickCheck 9 Test Cup  
 4215-DUA-187A-125 QuickCheck 8 Test Cup  
 4215-DUA-187A-150 QuickCheck 8 Test Cup  
 4215-DUA-177A-125 QuickCheck 7 Test Cup  
 4215-DUA-177A-150 QuickCheck 7 Test Cup