BioStrip HCG

One Step Urine Pregnancy test

REF 4220-5

INTENDED USE

For the rapid determination of human chorionic gonadotropin (hCG) in urine specimens. This test kit is used to obtain a visual, qualitative result for early detection of pregnancy. This product is for professional use only.

SUMMARY AND EXPLANATION

Human chorionic gonadotropin (hCG) is a glycoprotein hormone secreted by the developing placenta shortly after implantation. hCG can be detected in the urine and serum of pregnant women as early as 6 to 15 days after conception. The concentration of hCG increases to 5-50 IU/L one week post implantation, reaches about 100 IU/L at the time of the first missed menstrual period, and peaks at 100,000-200,000 IU/L at the end of the first trimester¹.

The appearance of hCG soon after conception and its subsequent rise in concentration during early gestational growth make it an excellent marker for the early detection of pregnancy.

The BioStrip HCG test kit is a rapid test to detect the presence of hCG in urine specimens in a qualitative format sensitive to 20 IU hCG/L. The test utilizes a combination of monoclonal and polyclonal antibody reagents to selectively detect elevated levels of hCG in urine. The immunological specificity of the test kit virtually eliminates cross reactivity interferences from the structurally related glycoprotein hormones hFSH, hLH, hTSH at physiological levels.

PRINCIPLE

The BioStrip HCG test is a qualitative, two site sandwich immunoassay for the determination of hCG in urine. The membrane was pre-coated with goat anti-hCG on the test band region and goat anti-mouse on the control band region. During the test, the patient urine is allowed to react with the coloured colloidal gold-mono-clonal anti-hCG conjugate which was pre-dried on the test strip. The mixture then moves upward on the membrane by the capillary action. For a positive result, a colour band with the specific antibody-hCG-coloured conjugate complex will form on the membrane. Absence of this coloured band in the test band region suggest a negative result. To serve as a procedural control, a coloured band at control region will always appear regardless the presence of hCG in test specimen.

SPECIMEN COLLECTION

Any urine specimen is appropriate for hCG testing, but the first morning urine is optimal because it generally contains the highest concentration of hCG. Urine specimens may be collected in any clean, dry, plastic or glass container. Specimens may be kept at room temperature for 8 hours or refrigerated at 2 - 9 °C for up to 72 hours. Samples may be frozen at -20 °C or below. If frozen, mix after thawing. Do not refreeze. Do not add preservatives. Specimen should be allowed to reach room temperature before testing. Very turbid urine specimens should be centrifuged or filtered prior to use. Samples with visible precipitates should be allowed to settle and the clear supernatant used for testing.

REAGENTS AND MATERIALS SUPPLIED

1. Test strips with goat anti-hCG and coloured colloidal gold – monoclonal mouse anti-hCG conjugate.
2. Package insert

MATERIALS REQUIRED BUT NOT SUPPLIED

1. Specimen collection container
2. Timer

PRECAUTIONS

1. For in vitro diagnostic use only.
2. Do not use after the expiration date
3. The cap of the container houses a desiccant to keep the test strips dry. After removing each strip, immediately replace the cap and close tightly to protect from light and moisture.
4. Store strips in the closed can until use. Perform testing immediately after removing each test strip from the can.
5. All strips must remain in the original can. Transfer to another container may cause test strips to deteriorate and become non-reactive.
6. Do not touch the test area of the strip.
7. Urine specimens and used strips should be considered potentially hazardous. Handle in the same manner as an infectious agent. Dispose according to local regulations.
8. For single use only. Do not reuse.


STORAGE AND STABILITY

Store at a temperature between 2°C to 30°C in the desiccated container with the lid firmly closed. Keep away from heat, moisture, and direct sunlight. Test strips are stable until the expiry date printed on the can.

TEST PROCEDURE

Bring test strips, patient samples, controls and reference materials to room temperature 20°C to 30°C prior to testing.


1. Remove the test strip from the desiccated container. Cap the container tightly and promptly immediately after use. Label the strip with patient or control identification.
2. Immerse the strip into the urine with the arrow end pointing towards the urine. Do not immerse past the MAX (maximum) line. You may leave the strip in the urine or you may take the strip out after a minimum of 10 seconds and lay the strip flat on a clean, dry, non-absorbent surface (e.g. mouth of the urine container).
3. Coloured band(s) will appear. Depending on the concentration of hCG in the test specimen, positive results may be observed as soon as 10-30 seconds. However, confirm negative results, the complete reaction time of 5 minutes is required. **Do not interpret results after 10 minutes.**

INTERPRETATION OF RESULTS

1. **Positive:** Distinct colour bands appear on the control and test region. Colour intensity of the test bands may vary.
2. **Negative:** Only one colour band appears on the control region. No visible band on the test region.

3. **Invalid:** No visible band at all or no visible band at the control region. Repeat test.

QUALITY CONTROL

In internal procedural control is included in the test. The colour band appearing on the control region of the strip indicates proper performance and reactive reagents.

Good laboratory practice suggest the use of external controls to ensure proper procedure and performance. Negative and positive controls containing hCG at various concentrations are available commercially.

LIMITATIONS

1. A number of conditions other than pregnancy, including trophoblastic disease and certain nontrophoblastic neoplasms, cause elevated levels of hCG. These diagnosis should be considered if appropriate to the clinical evidence.
2. If the urine specimen is too dilute (i.e. low specific gravity) it may not contain representative level of hCG. If pregnancy is still suspected, a first morning urine should be obtained from the patient 48-72 hours later and tested.
3. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after clinical and laboratory findings have been evaluated.
4. False negative reaction may occur when the HCG level in the urine specimen is extremely high. This typically occurs when the patient's HCG levels peaks at 13-16 week pregnancy. If this is suspected, one can dilute one drop of urine with 50 drops of saline and retest the diluted urine using a fresh test strip. If the diluted specimen gives a positive result, the pregnancy test can be read as positive.

PERFORMANCE CHARACTERISTICS

Specificity

The specificity of the BioStrip HCG test device was evaluated in an experiment to determine the cross reactivity of structurally related hormones: human luteinizing hormone (hLH), human follicle stimulating hormone (hFSH), and human thyroid stimulating hormone (hTSH).

	hLH 500 IU/L	hFSH 1000 IU/L	hTSH 1000 mIU/L
Positive	0/20	0/20	0/20
Negative	20/20	20/20	20/20

BioStrip HCG testing urine samples spiked with 0 IU/L hCG and structurally related hormones

	hLH 500 IU/L	hFSH 1000 IU/L	hTSH 1000 mIU/L
Positive	20/20	20/20	20/20
Negative	0/20	0/20	0/20

BioStrip HCG testing urine samples spiked with 20 IU/L hCG and structurally related hormones

Sensitivity

The analytical sensitivity of BioStrip HCG test was validated by an experiment to determine the optimal urine hCG concentrations. All 0 IU/L specimens showed negative results. At 10 IU/L, 24 of 40 specimens tested positive. At 15 IU/L, 36 of 40 specimens tested positive. All samples ≥ 20 IU/L showed positive results within 5 minutes. This study indicates that BioStrip HCG test can consistently detect 20IU/L hCG in urine specimens, showing that 20IU/L is the detectable sensitivity for BioStrip tests. Results are tabulated in the table below.

hCG IU/L	0	10	15	20	25	50	100
Positive	0	24	36	40	40	40	40
Negative	40	16	4	0	0	0	0
Total n	40	40	40	40	40	40	40

BioStrip HCG testing urine samples spiked with various concentrations of hCG, where n=40

Standardization

The BioStrip HCG test has been standardized to World Health Organization Third International Standard (WHO 3rd IS 75/589).

Interference Testing

The following substances were added to hCG-free urine samples and urine samples spiked with 25 IU hCG/L. None of the substances at the listed concentrations interfered with the assay results.

Acetaminophen	20 mg/mL
Acetylsalicylic acid	20 mg/mL
Ascorbic acid	20 mg/mL
Atropine	20 mg/mL
Caffeine	20 mg/mL
Gentesic acid	20 mg/mL
Glucose	2 g/dL
Haemoglobin	1 mg/dL

REFERENCES

1. Batzer FR, Fertility & Sterility, Vol. 34, 1, 1980.
2. Catt KJ, Dufan ML & Vaitukaitis JL, J. Clin. Endocrinol. Metab., Vol 40, 537, 1975.
3. Braunstein GD, Rasor J, Alder D, Danzer H, Wade ME, Am.J. Obstet. Gynecol., Vol 126, 678, 1976.
4. Lenton EA, Neal LM, Sulaiman R, Fertility & Sterility, Vol. 37, 773, 1982.
5. Dawood MY, Saeba BB, & Landesman R, Ob. Gyn. Vol. 126, 678, 1976.
6. Braunstein GD, et al. Ann. Inter. Med. Vol. 78, pp 39-45, 1973.
7. Engvall E, Methods in Enzymology, Vol. 70, pp 419-439, 1980.
8. Uotila M, Ruoslahti E, Engvall EJ, Immunol. Methods, Vol. 42, 11, 1981

Advena Ltd
Pure Offices, Plato Close
Warwick CV34 6WE UK

Innovatek Medical Inc.
3 - 1600 Derwent Way, Delta BC Canada V3M 6M5

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