



REF COV-19C25

COVID-19 Antigen Rapid Test Device

Health Canada Approved - Point of Care

Specimen: Nasal / Nasopharyngeal Swab




Our new **Rapid Response® COVID-19 Antigen Rapid Test Device** allows for an efficient and cost-friendly method of detection. In just 15 minutes, the test can detect the SARS-CoV-2 viral nucleoprotein antigen through visual interpretation of colour development. The test is intended for use with Nasal and Nasopharyngeal swab specimens within 6 days of symptom onset.

Easy to Use:

1. Gently mix buffer. Add 10 drops of buffer to the tube.
2. Collect nasal/nasopharyngeal specimen with the swab and then insert swab into the tube. Swirl the swab and mix well.
3. Squeeze the swab against the inside of the tube 10-15 times to release all the liquid. Dispose of the swab properly.
4. Attach nozzle onto the tube.
5. Add 3 drops to the sample well.
6. **See results in 15 minutes!**

Clinical Evaluation - Nasal Swab (Anterior Nares)




Diagnostic Sensitivity
CI: 78.6% - 96.2%



Diagnostic Specificity
CI: 96.5% - 100%

Clinical Evaluation - Nasopharyngeal Swab




Diagnostic Sensitivity
CI: 79.8% - 99.3%



Diagnostic Specificity
CI: 84.5% - 100%

Interpretation of Test Results



Positive



Negative



Invalid

Kit Contents

- Individually Packed Test Devices
- Individually Packed Swabs
- Tube Stand
- Extraction Tube
- Extraction Buffer
- Nozzle with Filter
- Package Insert



Tel: +1 905-944-9565
Email: covid19@btnx.com
www.btnx.com

Website
btnx.com



Product Insert
rb.gy/7boeas




Molecular VS Rapid Response® COVID-19 Antigen Rapid Test Device

	Molecular Testing (RT-PCR)	Rapid Response® COVID-19 Antigen Test
Principle	Nucleic Acid Test of COVID-19	SARS-CoV-2 Viral Nucleoprotein Antigen
Time to Results	> 6 hours	15 minutes
Specimen	Throat, Nasopharyngeal Fluid	Nasal, Nasopharyngeal Swab
Test Capacity	Limited Throughput	High Turnover
Characteristics	<ul style="list-style-type: none"> - Detects active infection in patients - Testing requires equipment primarily housed in laboratory settings - Does not identify prior infections 	<ul style="list-style-type: none"> - Nasal / Nasopharyngeal Swab: Detects within 6 days of symptom onset - Does not require specialized laboratory equipment - Rapid results can be easily interpreted

SARS-CoV-2 Testing Timeline

The chart below gives a representation of SARS-CoV-2 Infection and the onset of antibodies in an infected individual.



Ordering Information

Product Code	Product Name	Contents
F27-COV-19C25	COVID-19 Antigen Rapid Test Device	25 Tests / Kit