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Panel Height
11.000

LABELING ARTWORK DATA
BAXTER CONFIDENTIAL - INTERNAL USE ONLY
Baxter Healthcare Corporation , Medical Products 17511 Armstrong Avenue, Irvine, California 92614-5725 USA Telephone (949) 851-9066, Fax (949) 474-6455
Labeling Artwork Filename: Labeling Artwork 071971058.indd
Labeling Artwork Electronic File Format: Adobe InDesign CS6
Labeling Artwork Page Size: 11" wide x 17" high
Labeling Artwork Revised By: Dan Futoran Labeling Artwork Revision Date: 16-Jul-2013
Unless otherwise specified, all dimensions are shown as fractional inches with decimal equivalents in parentheses.

LABELING ARTWORK PROOFREADING APPROVALS	Language: French
Electronic Proofreading Approvals <input type="checkbox"/> N/A (section shall remain blank when checked), Initials/Date:	
<input type="checkbox"/> See attached	
Inspection Report Printed Name _____ Signature _____ Date _____	
Manual Proofreading Approvals <input type="checkbox"/> N/A (section shall remain blank when checked), Initials/Date:	
<input type="checkbox"/> 1st Proofread Printed Name _____ Signature _____ Date _____	
<input type="checkbox"/> 2nd Proofread Printed Name _____ Signature _____ Date _____	
IRF03004, REV C	


Language Legend	
ENG	English
FRE	French
GER	German
SPA	Spanish
ITA	Italian
DUT	Dutch
SWE	Swedish
NOR	Norwegian
DAN	Danish
FIN	Finnish
POR	Portuguese
GRE	Greek

Title: **DIRECTION INSERT, CE INTERMATE SYSTEM, BASIC**

Panel 13 — Flow Profile panel

Panel Width **8.500**

INTERMATE

1 ENG) Incremental Flow Profile
 FRE) Profil incrémentiel du débit
 GER) Flussprofil
 SPA) Incremento del perfil de flujo
 ITA) Profilo incrementale di flusso
 DUT) Inloopprofiel met periodieke verhogingen
 SWE) Stigande flödesprofil
 NOR) Profil for stigende infusjonshastighet
 DAN) Profil for stigende infusionshastighed
 FIN) Nouseva virtausprofiili
 POR) Perfil de Débito Incremental
 GRE) Προφίλ αυξανόμενης ροής

2 ENG) Percent of Nominal Flow Rate (ml/h)
 FRE) Pourcentage du débit nominal (ml/h)
 GER) Prozent der nominalen Flussrate (ml/h)
 SPA) Porcentaje de velocidad de flujo nominal (ml/h)
 ITA) Percentuale della velocità di flusso nominale (ml/h)
 DUT) Percentage van nominale inloopsnelheid (ml/uur)
 SWE) Procent av nominell flödes hastighet (ml/h)
 NOR) Prosent av nominell infusjonshastighet (ml/t)
 DAN) Procent af nominel flowhastighed (ml/t)
 FIN) Prosenttia nimellisestä virtausnopeudesta (ml/h)
 POR) Percentagem de Taxa de Débito Nominal (ml/h)
 GRE) Ποσοστό ονομαστικού ρυθμού ροής (ml/h)

3 ENG) Percent of Nominal Delivery Time (hours)
 FRE) Pourcentage du temps d'administration nominal (heures)
 GER) Prozent der nominalen Verabreichungsdauer (Stunden)
 SPA) Porcentaje de tiempo de perfusión nominal (horas)
 ITA) Percentuale del tempo di erogazione nominale (ore)
 DUT) Percentage van nominale toedieningstijd (uur)
 SWE) Procent av nominell administreringstid (timmar)
 NOR) Prosent av nominell administreringstid (timer)
 DAN) Procent af nominel infusionstid (timer)
 FIN) Prosenttia nimellisestä virtausajasta (tunteja)
 POR) Percentagem de Tempo de Administração Nominal (h)
 GRE) Ποσοστό ονομαστικού χρόνου χορήγησης (ώρες)

Panel Height **11.000**

LABELING ARTWORK DATA
BAXTER CONFIDENTIAL - INTERNAL USE ONLY
 Baxter Healthcare Corporation, Medical Products
 17511 Armstrong Avenue, Irvine, California 92614-5725 USA
 Telephone (949) 851-9066, Fax (949) 474-6455
 Labeling Artwork Filename: **Labeling Artwork 071971058.indd**
 Labeling Artwork Electronic File Format: Adobe InDesign CS6
 Labeling Artwork Page Size: 11" wide x 17" high
 Labeling Artwork Revised By: **Dan Futoran** Labeling Artwork Revision Date: **16-Jul-2013**
Unless otherwise specified, all dimensions are shown as fractional inches with decimal equivalents in parentheses.

LABELING ARTWORK PROOFREADING APPROVALS		Language: Multilingual—12 languages	
Electronic Proofreading Approvals <input type="checkbox"/> N/A (section shall remain blank when checked), Initials/Date:			
<input type="checkbox"/> See attached	Printed Name	Signature	Date
Manual Proofreading Approvals <input type="checkbox"/> N/A (section shall remain blank when checked), Initials/Date:			
<input type="checkbox"/> 1st Proofread	Printed Name	Signature	Date
<input type="checkbox"/> 2nd Proofread	Printed Name	Signature	Date

Language Legend	
ENG	English
FRE	French
GER	German
SPA	Spanish
ITA	Italian
DUT	Dutch
SWE	Swedish
NOR	Norwegian
DAN	Danish
FIN	Finnish
POR	Portuguese
GRE	Greek

Title: **DIRECTION INSERT, CE INTERMATE SYSTEM, BASIC**

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

Panel Height **11.000**

LABELING ARTWORK DATA
BAXTER CONFIDENTIAL - INTERNAL USE ONLY
 Baxter Healthcare Corporation, Medical Products
 17511 Armstrong Avenue, Irvine, California 92614-5725 USA
 Telephone (949) 851-9066, Fax (949) 474-6455
 Labeling Artwork Filename: **Labeling Artwork 071971058.indd**
 Labeling Artwork Electronic File Format: **Adobe InDesign CS6**
 Labeling Artwork Page Size: **11" wide x 17" high**
 Labeling Artwork Revised By: **Dan Futoran** Labeling Artwork Revision Date: **16-Jul-2013**
Unless otherwise specified, all dimensions are shown as fractional inches with decimal equivalents in parentheses.

Language: **Multilingual—12 languages**

LABELING ARTWORK PROOFREADING APPROVALS			
Electronic Proofreading Approvals <input type="checkbox"/> N/A (section shall remain blank when checked), Initials/Date:			
<input type="checkbox"/> See attached	_____	_____	_____
Inspection Report	Printed Name	Signature	Date
Manual Proofreading Approvals <input type="checkbox"/> N/A (section shall remain blank when checked), Initials/Date:			
<input type="checkbox"/> 1st Proofread	_____	_____	_____
	Printed Name	Signature	Date
<input type="checkbox"/> 2nd Proofread	_____	_____	_____
	Printed Name	Signature	Date

IRF03004, REV B


Language Legend	
ENG	English
FRE	French
GER	German
SPA	Spanish
ITA	Italian
DUT	Dutch
SWE	Swedish
NOR	Norwegian
DAN	Danish
FIN	Finnish
POR	Portuguese
GRE	Greek

Title: **DIRECTION INSERT, CE INTERMATE SYSTEM, BASIC**

Panel 15 — INTERMATE SV and INTERMATE LV illustration panel


Panel Width **8.500**

INTERMATE SV

Panel Height **11.000**

INTERMATE LV

LABELING ARTWORK DATA
BAXTER CONFIDENTIAL - INTERNAL USE ONLY
 Baxter Healthcare Corporation, Medical Products
 17511 Armstrong Avenue, Irvine, California 92614-5725 USA
 Telephone (949) 851-9066, Fax (949) 474-6455
 Labeling Artwork Filename: **Labeling Artwork 071971058.indd**
 Labeling Artwork Electronic File Format: **Adobe InDesign CS6**
 Labeling Artwork Page Size: **11" wide x 17" high**
 Labeling Artwork Revised By: **Dan Futoran** Labeling Artwork Revision Date: **16-Jul-2013**
Unless otherwise specified, all dimensions are shown as fractional inches with decimal equivalents in parentheses.

LABELING ARTWORK PROOFREADING APPROVALS		Language: Multilingual—12 languages	
Electronic Proofreading Approvals <input type="checkbox"/> N/A (section shall remain blank when checked), Initials/Date:			
<input type="checkbox"/> See attached			
Inspection Report	Printed Name	Signature	Date
Manual Proofreading Approvals <input type="checkbox"/> N/A (section shall remain blank when checked), Initials/Date:			
<input type="checkbox"/> 1st Proofread	Printed Name	Signature	Date
<input type="checkbox"/> 2nd Proofread	Printed Name	Signature	Date


Language Legend	
ENG	English
FRE	French
GER	German
SPA	Spanish
ITA	Italian
DUT	Dutch
SWE	Swedish
NOR	Norwegian
DAN	Danish
FIN	Finnish
POR	Portuguese
GRE	Greek

Title: **DIRECTION INSERT, CE INTERMATE SYSTEM, BASIC**

Panel 16 — Back panel/INTERMATE XLV illustration panel


Panel Width **8.500**

INTERMATE XLV

Panel Height **11.000**

Baxter Healthcare SA
 8010 Zurich
 Switzerland

Baxter Healthcare Corporation
 Deerfield IL 60015 USA

Made in USA

* + H 1 6 0 0 7 1 9 7 1 0 5 8 2 J *

LABELING ARTWORK DATA

BAXTER CONFIDENTIAL - INTERNAL USE ONLY

Baxter Healthcare Corporation, Medical Products
 17511 Armstrong Avenue, Irvine, California 92614-5725 USA
 Telephone (949) 851-9066, Fax (949) 474-6455

Labeling Artwork Filename: **Labeling Artwork 071971058.indd**

Labeling Artwork Electronic File Format: **Adobe InDesign CS6**

Labeling Artwork Page Size: **11" wide x 17" high**

Labeling Artwork Revised By: **Dan Futoran** Labeling Artwork Revision Date: **16-Jul-2013**

Unless otherwise specified, all dimensions are shown as fractional inches with decimal equivalents in parentheses.

LABELING ARTWORK PROOFREADING APPROVALS

Language: **Multilingual—12 languages**

Electronic Proofreading Approvals N/A (section shall remain blank when checked), Initials/Date:

See attached

Inspection Report Printed Name _____ Signature _____ Date _____

Manual Proofreading Approvals N/A (section shall remain blank when checked), Initials/Date:

1st Proofread

Printed Name _____ Signature _____ Date _____

2nd Proofread

Printed Name _____ Signature _____ Date _____