


SPEEDSPLINT™

MULTI-PURPOSE SPLINT



*Fast & Simple
Temporary
Immobilization
Low Cost*



FEMORAL TRACTION KIT

*Easily Applied in
3 to 4 Minutes
Lightweight
High-Tech*



SPEEDSPLINT™

MULTI-PURPOSE SPLINT

RE-ORDER #M100

SIMPLE TEMPORARY IMMOBILIZATION ANYONE CAN APPLY WITHOUT THE RISK OF SHUTTING OFF THE BLOOD SUPPLY

● EMERGENCY CARE PROFESSIONALS

Supports good padding procedures & scissor trims & tapes for customized patient fit.

Lightweight compact flat for efficient use of limited & valuable space.

● X-RAYS PENETRATE

● ONE SIZE FITS ALL


Very versatile, can be used on foot, ankle, lower leg, knee, hand, wrist, lower arm, elbow or upper arm.

● VOLUNTEER AMATEUR SPORTS COACHES

Fast, simple application on the playing field, restricts movement & protects the player when needed.

● CHILDREN'S FRACTURES & SPRAINS

About 4 out of 5 cases transport to the hospital by family, neighbour car or taxi with mother, in first response care....often more unnecessary injury at this time because of either inadequate or no simple immobilization.



CONVENIENT

TO SHIP • TO STORE • TO REPLACE

SPEEDSPLINT™

FEMORAL TRACTION KIT

RE-ORDER #F200

KIT CONTAINS: • FEMORAL TRACTION SPLINT • PADDING • TRACTION HARNESS • GROIN PADS

● SPACE STINGY

● PRECISE TRACTION TENSION CONTROL

● X-RAYS PENETRATE

Fits any size thigh and up to a 6'4" or 1.73m person

Just 2.2 lbs. or 1 Kg., this complete emergency response plastic sealed carton efficiently saves precious disaster storage and rescue vehicle space . . . only 20" x 15" x 2-1/2" or 51 cm x 38 cm x 6.5 cm . . . floats too!

Speedsplint Splints are manufactured from an exceptionally high-tech, light-weight, strong single piece EXTRUDED "yellow alert" plastic . . . and it's Re-usable!

2 EMT's can secure "fingertouch" Velcro traction control, exerting NO proximal shaft flexion, maintaining easy access to monitor and adjust traction tension during transport.

After traction is secured by emergency response team, only the Orthopedic Surgeon need remove the patient from traction

™ SPEEDSPLINT IS A REGISTERED TRADEMARK

MANUFACTURED BY:

GEEN
HEALTH CARE
I N C .

931 PROGRESS AVENUE • UNIT 13
SCARBOROUGH • ONTARIO
CANADA • M1G 3V5
TOLL FREE: 1-800-565-GEEN
TELEPHONE: 416-439-2237
FAX: 416-439-6934

DISTRIBUTED BY:

