



Non-Adherent Dressings







Pain-free removal.

Telfa Non-Adherent Dressings

- Feature a perforated polyester film, containing a dense, non-woven cotton absorbent core that can be used on a variety of wounds
- Have a non-adherent surface with smooth aperture edges for intact and pain-free removal
- Allow exudate to pass through optimal size apertures without allowing tissue to adhere to the dressing
- Will not adhere to sutures, staples or adhesive wound closures


Under the Microscope

Here you can see how Telfa dressing's superior apertures compare to our competitors.¹ Notice the optimal size apertures of the Telfa dressing, which allow exudate to pass through without allowing tissue to adhere to the dressing.



A Comparison of Adherency Among Dry, Non-Adherent Dressings

In this test, adherence was measured by the amount of force required to remove the dressing from a two-inch strip of Bovine liver.¹




Source:-Quigg, J. Gilbert, J & Reil, L 1994



Telfa Clear Dressings

Non-Adherent perforated polyester film


Non-Adherent	Allows for new tissue growth without interrupting the healing process. Allows for pain-free removal; Will not adhere to sutures, staples, or adhesive wound closures
Non-Absorbent	When combined with an appropriate secondary dressing, may be used on lightly, moderately, or heavily draining wounds
Clear Film	Allows assessment of the wound site without disrupting the healing process; May be left in place for several days
Versatile	Available in many convenient sizes; May be used on a variety of wounds
May be Cut	Allows for custom wound applications



Telfa "Ouchless" Non-Adherent Dressings

Non-adherent perforated polyester film, plus absorbent pad


Non-Adherent	Allows for new tissue growth without interrupting the healing process; Allows for pain-free removal; Will not adhere to sutures, staples, or adhesive wound closures
Absorbent Core	Cotton core offers light to moderate absorbency to remove excess fluid from wound site
Versatile	Available in many convenient sizes; May be used on a variety of wounds



Telfa Island Dressings and Telfa Plus Island Dressings

Non-Adherent perforated polyester film, plus pad with adhesive border

Non-Adherent	Allows for new tissue growth without interrupting the healing process; Allows for pain-free removal; Will not adhere to sutures, staples, or adhesive wound closures
Absorbent Core	Cotton core offers light to moderate absorbency to remove excess fluid from wound site (Telfa Plus offers thicker pad for greater absorbency)
Strike-Through Barrier	Unique adhesive border helps prevent fluid strike-through and resists external bacteria penetration
Soft Non-woven Backing	Conformable; comfortable to wear
Adhesive Border	No requirement for secondary dressing for fixation; Saves nurses' time. Ideal dressing for surgical wounds
Versatile	Available in many convenient sizes; May be used on a variety of wounds



Telfa Clear Wound Dressings

Code	Description	Per Box	Ship Case
Sterile 1's in Peel-Back Package			
1109	3" x 3"	25	50
1111	4" x 5"	25	50
1113	12" x 12"	25	50
1114	12" x 24"	25	50

Telfa Clear Wound Dressings (Bulk)

Code	Description	Ship Case
Nonsterile Bulk		
1115	39" x 25 yds.	4 rolls

Telfa "Ouchless" Non-Adherent Dressings

Code	Description	Per Box	Ship Case
Sterile 1's in Peel-Back Package			
1050	3" x 4"	50	900
1169	3" x 6"	50	750
1238	3" x 8"	50	600
28959	8" x 10"	25	400

Sterile 1's in Peel-Back Package, Perforated Carton for Easy Dispensing

1961	2" x 3"	100	2,400
2132	3" x 4"	100	2,400

Nonsterile Bulk, Non-Adherent Film on One Side

2891	3" x 8"	200	1,000
3279	8" x 10"	125	500

Nonsterile Bulk, Non-Adherent Film on Both Sides

12388	3" x 8"		4,000
-------	---------	--	-------

Telfa Sterile Strips

Code	Description	Per Box	Ship Case
Sterile Telfa Dressing			
28959	8" x 10"	25	400

Telfa Island Dressings

Code	Description	Per Box	Ship Case
Sterile 1's in Peel-Back Package			
7539LF	2" x 3¾"	50	400
7540	4" x 5"	25	200
7541	4" x 8"	25	100
7542	4" x 10"	25	100
7544	4" x 14"	25	50
7550	4" x 4"	25	200
7551	6" x 6"	25	100

Telfa Plus Island Dressings

Code	Description	Per Box	Ship Case
Sterile 1's in Peel-Back Package			
2562	4" x 6"	25	100
2563	6" x 7"	25	100
2564	6" x 10"	25	100
2565	8" x 8"	15	60

Telfa "Ouchless" Adhesive Dressings

Code	Description	Per Box	Ship Case
Sterile 1's in Peel-Back Package			
6017	2" x 3"	100	2,400
7643	3" x 4"	100	1,200

¹ Quigg, J et al, Technical and Clinical Experience with Telfa Clear, 1994, Kendall Healthcare Products Co.

COVIDIEN, COVIDIEN with logo, Covidien logo and positive results for life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. © 2011 Covidien.

Manufactured by: Covidien Inc, 15 Hampshire Street, Mansfield, MA, USA, 02048

CDN-P00175-E (H6313) Rev. 2011/11

