# **ResMed**





Quick to fit. Effective to use.







#### PRODUCT CODE

ResMed Hospital Full Face Mask (vented)		ResMed Hospital Nasal Mas
Small	60703	Medium
Medium	60704	Large
Large	60705	Headgear (10 per pack)
Headgear (10 per pack)	60710	



61118

21952
21953

ResMed Corp Poway, CA, USA + 1858 746 2400 OR 1 800 424 0737 (toll free), ResMed Ltd. 1 Elizabeth MacArthur Drive Bella Vista NSW 2153 Australia +61 (2) 8884 1000 or 1 800 658 189 (toll free), Offices in Austria, Brazil, Finland, France, Germany, Hong Kong, Japan, Malaysia, Netherlands, New Zealand, Singapone, Spain, Sweden, Switzerland, United Kingdom (see website for details), Hospital Full Face Mask potentical by patents: AU 744539, AU 779557, AU 75975, CA 2573 (PD) C 9056069, FT 10756, US 551754, ResMed, a leading global supplier of masks, has developed a line of custom-designed products for hospital use. Quick to fit and easy to use, ResMed's disposable masks and circuits provide cost-effective, time-saving solutions for the clinical environment. Ideally suited for CPAP, bilevel pressure treatment and noninvasive ventilation therapy.


### **ResMed Hospital Full Face Mask**

The proven comfort and seal of ResMed's Mirage" cushion technology in a hospital mask. Quick to fit and easy to adjust, the full face mask is available in three sizes—**small, medium** and **large**—providing a comfortable fit for a wide range of patients. Designed for single patient use, the masks can be easily discarded, eliminating cross-contamination concerns and costs associated with mask disinfection and sterilization.



## **ResMed Hospital Nasal Mask**

ResMed's nasal masks are designed for single patient use. They are available in two sizes—**medium** and **large** providing a comfortable fit for a wide range of patients. The masks can be easily discarded, eliminating crosscontamination concerns and costs associated with mask disinfection and sterilization.



### **ResMed Universal Disposable Bilevel Circuits**

The universal bilevel circuits are designed for single patient use only and may be used with Res/Med Hospital full face and nasal masks. The circuits are available both with a filter (21952) and without a filter (21953). Each circuit is packaged with three additional components. Please refer to the user's guide for operating instructions.





72" Disposable Bilevel Circuit with filter (21952)

72" Disposable Bilevel Circuit (21953)

Additional components packaged with Bilevel Circuits



96" pressure monitoring tubing

30" "dual" pressure monitoring extension tubing and T-assembly Inline exhalation port connector

Pressure sensing

### Saving you money... saving you time

#### So easy to use, so quick to fit

- Preassembled and ready to use straight from the package
- Color-coded labeling for easy size identification
- Fitting Guide on packaging ensures correct size selection
- Robust one-piece design—no assembly required
- Simple and thoughtful design with minimal adjustment points

#### Comfortable and effective therapy

- Effective CO<sub>2</sub> washout
- Reliable therapy custom-designed for the clinical environment
- ResMed sealing technology provides optimal comfort, leak-free seal and reduction of
- pressure sores
  A broad range of cushion sizes offers a comfortable fit

