

**EntraFlo™ Nutrition Delivery Systems**Providing Accurate, Reliable Tube Feeding

## Your Patient Requires Tube Feeding. Are You Well-Equipped to Meet Their Needs?

Tube feeding isn't something you encounter every day. But when your patient needs it, the right equipment can make all the difference. Because the last thing you want is to figure out how to use a complicated piece of equipment. Or troubleshoot one that's prone to breaking down.

We understand. And we can help.

#### Rely on a System That's Tried, Tested and Proven.

Our EntraFlo™ nutrition delivery systems make tube feeding easier for you and your patients.

- » Simple to set up and use
- » Accurate formula delivery
- » Intuitive operation

Formerly known as COMPAT® and COMPAT DualFlo®, these pumps have been the standard for enteral feeding for over 20 years.

## New Name. Same Solid Performance.



# Complete and Completely Reliable.

When your nutrition delivery system is user friendly and dependable, everyone benefits. Nurses are more productive. Patients receive better care. And families are more confident in your facility. Our EntraFlo pumps and administration sets have features that make operation and formula administration easy.

# Which Feeding System Is Right For You? With a few quick questions, our Feeding Pump & Set Selector can help you find the device that best meets your needs. To get started, visit medline.com/entraflo.







#### **EntraFlo Nutrition Delivery System**

This trusted tube feeding pump is simple to set up and adjust.

- » Dual display windows let nurses monitor infusion rate, volume delivered and dose limit simultaneously
- » Automatically stops when dose is finished
- » Remembers flow rate, dose limit, volume delivered and accumulated volume
- » 8-hour battery life at 100 ml/hr, rechargeable while pump is operating

Description

» Simple to set up\*

Itam No

#### EntraFloH<sub>2</sub>O™ Nutrition Delivery System

This advanced delivery system can be programmed for individualized patient feeding and hydration.

- » Auto primes formula and water through the administration set
- » Security doors restrict access to the administration set tubing and help keep it from being knocked loose
- » Remembers flow rate, dose limit, formula strength, water volume set, time interval, formula and water volume delivered, and accumulated volume
- » 8-hour battery life at 100 ml/hr; rechargeable while pump is operating
- » Simple to set up\*

Dla


item No.	Description	PKg.
ENT199235RNT	EntraFlo Feeding Pump, Rental	1/ea
ENT199235L	EntraFlo Feeding Pump, Lease	1/ea
ENT199235	EntraFlo Feeding Pump, Purchase	1/ea
ENFIT70550	EntraFlo Safety Spike Plus Set, ENFit Connector	30/cs
ENFIT70552	EntraFlo 1000 ml Bag Set, ENFit Connector	30/cs
ENFIT1000TC	ENFit Transition Connector for ENFit Syringe or ENFit Feeding Set, Connects to ENFit Feeding Tube	30/cs
ENFIT1010GC	ENFit Transition Connector for Non- ENFit Syringe or Non-ENFit Feeding Set, Connects to Non-ENFit Feeding Tube	30/cs

Item No.	Description	Pkg.
ENT199255RNT	EntraFloH₂O Feeding Pump, Rental	1/ea
ENT199255L	EntraFloH₂O Feeding Pump, Lease	1/ea
ENT199255	EntraFloH₂O Feeding Pump, Purchase	1/ea
ENFIT70554	EntraFloH₂O Safety 1000 ml Feed and 1000 ml Flush Bag Set, ENFit Connector	30/cs
ENFIT70555	EntraFloH <sub>2</sub> O Safety Spike Plus and 1000 ml Flush Bag Set, ENFit Connector	30/cs
ENFIT1000TC	ENFit Transition Connector for ENFit Syringe or ENFit Feeding Set, Connects to ENFit Feeding Tube	30/cs
ENFIT1010GC	ENFit Transition Connector for Non- ENFit Syringe or Non-ENFit Feeding Set, Connects to Non-ENFit Feeding Tube	30/cs

<sup>\*</sup>Illustrations and complete instructions for use available at medline.com/entraflo

### **Delivering Optimal Performance and Value.**

You always put your patients first. So do we. We offer a variety of gastrostomy tubes and feeding syringes to help you provide comfortable and safe enteral care.



#### **Medline Gastrostomy Tubes**

Rounded tips ensure that our gastrostomy tubes insert easily. Closure caps seal the tube between feedings.

- » Gastric balloon expands evenly for a secure, comfortable fit
- » Circular skin disk easily adjusts to maintain proper tube position, plus allows stoma observation
- » Made of 100 percent medical-grade silicone

Item No.	Description	Balloon Size	Pkg.
DYND70316	White, 16 Fr	15 cc	1/ea
DYND70318	White, 18 Fr	15 cc	1/ea
DYND70320	White, 20 Fr	15 cc	1/ea
DYND70322	White, 22 Fr	15 cc	1/ea
DYND70324	White, 24 Fr	15 cc	1/ea

#### **Medline Feeding Syringes**

Our 60 cc syringes have an elongated tip for smooth, accurate tube feeding and water flushes.

» Bold cc and ounce markings

Item No.	Description	Pkg.
DYND70642	Feeding Syringe, Includes Pole Bag and Luer Tip Adapter	30/cs
DYND20335	Feeding Syringe, Includes Pole Bag and 500 ml Graduated Container	30/cs
DYNC7061	Feeding Syringe, Includes Tray and 500 ml Graduated Container	30/cs
DYND70651	Bulb Syringe, Includes Pole Bag	30/cs
ENT6006M	Sterile Feeding Syringe	50/cs




#### **Making Safer Connections**

The ENFit connector reduces the risk of tube feeding misconnections, helping enhance patient safety. It has a unique enteral-specific design that:

- » Doesn't allow administration sets, gastrostomy tubes and syringes to connect with any other connector for any other clinical use
- » Has a locking feature that secures the connection

To learn more, visit medline.com/enfit.



#### **ENFit® Gastrostomy Tubes**

Made to comply with the new ENFit connection standard, these tubes feature an ergonomic cap for simple and secure closure.

- » Gastric balloon expands evenly for secure, comfortable fit
- » Circular skin disk easily adjusts to maintain proper tube position, plus allows stoma observation
- » Made from medical-grade silicone

Item No.	Description	Balloon Size	Pkg.
ENFIT70314	Clear, 14 Fr	5 cc	5/cs
ENFIT70316	Clear, 16 Fr	5 cc	5/cs
ENFIT70318	Clear, 18 Fr	20 cc	5/cs
ENFIT70320	Clear, 20 Fr	20 cc	5/cs
ENFIT70322	Clear, 22 Fr	20 cc	5/cs
ENFIT70324	Clear, 24 Fr	20 cc	5/cs

#### **ENFit Feeding Syringes**

These 60 cc syringes meet the ENFit connection standard and have a locking feature that secures the connection for accurate tube feeding.

» Bold cc and ounce markings

Item No.	Description	Pkg.
ENFIT70642	ENFit Feeding Syringe, Includes Pole Bag	30/cs
ENFIT20335	ENFit Feeding Syringe, Includes Pole Bag and 500 ml Graduated Container	30/cs
ENFIT7061	ENFit Feeding Syringe, Includes Tray and 500 ml Graduated Container	30/cs



COMMITTED AND READY TO SUPPORT YOUR NUTRITIONAL CARE.

Why choose Medline for your nutrition feeding systems? As the largest privately held manufacturer and distributor of med./surg. supplies in the U.S., we go above and beyond to give you exceptional quality, service and value. Always.

#### **Therapeutic Nutritional Support**

To help you meet the unique needs of your patients, we offer a wide variety of recognized and trusted nutrition formulas and modular supplements. The high protein content of Medline Active® Liquid Protein—15 to 21 grams for every 1-oz. serving—provides roughly one-quarter of the suggested daily value.

- » Ready-use-formula for tube feeding
- » Available in both sugar and sugar-free options

#### **Manufacturer-Direct Pricing**

Unlike other providers, we control the entire supply chain. We leverage economies of scale. Innovate operational performance. Eliminate third-party costs. And pass the savings straight to you.

- » Growing network of ISO-certified and quality-compliant manufacturing centers
- » Rigorous quality control and testing

#### **Hands-On Service**

Having a readily available Medline Representative and support network with the expertise to problem solve with you is priceless.

- » Proactive attention of a dedicated Medline Representative focused on your success
- » Customer service on call around the clock for overnight shipping and emergency orders
- » When you order a replacement pump, it'll arrive within 48 hours

#### **Ease of Conversion and Implementation**

Our experts will work with you to ensure your conversion to Medline EntraFlo is smooth and timely.

- » Clinical consultation and support
- » Product training webinars and in-services

#### **Speedy Delivery**

To make sure you get your supplies fast, our warehouses are located close to our customers.

- » 43 distribution centers in North America
- » 425+ trucks nationwide
- » 1-day ship to 96% of the U.S.

To find out more or to place an order, call your Medline Representative or Customer Service at 1-800-MEDLINE. In Canada, customers may call 1-800-396-6996. For additional resources, visit www.medline.com/entraflo.





Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

Medline United States 1-800-MEDLINE (633-5463) medline.com | info@medline.com Medline Canada 1-800-396-6996 medline.ca | canada@medline.com Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com