





# STRONG ADHESION. ACCURATE TRACINGS.

**MedGel® ECG Electrodes**Exceptional Reliability, Quality and Selection

# STRONG PERFORMANCE FOR ACCURATE TRACINGS EVERY TIME.

To best diagnose your patient, you need a clear ECG signal. But if the electrodes don't stay securely in place and don't have a high performing gel, you can't get it. We understand. And we can help. Our MedGel® Electrodes deliver exceptional adhesive strength and high performing gel for good conductivity. Have unique features that simplify work flow. And come in a full range of choices to meet your clinical priorities.

#### **Reliable Monitoring**

- » Strong adhesive holds fast to skin, yet removes cleanly
- » Highly conductive gel ensures optimal signal quality
- » Resealable bag helps keep electrodes from drying out (50 pack only)


#### Easy to Handle

- » Perforated strip backing provides quick, simple lead wire attachment
- » Full-edge lift tab allows easy removal
- » Gridlines on the strip backing help reduce clinician slips

#### **Comprehensive Selection**

- st Choose from a variety of materials, shapes and sizes
- » Color-coded packaging streamlines product selection





Gridlines on the strip backing promote clinician safety.



Resealable bag prolongs freshness (50 pack only).



Color-coded packaging makes product selection intuitive.

# **General Monitoring**

With their specially formulated high-performance gel and extended conductive surface, our general monitoring electrodes give you clear, reliable tracings.


#### **MedGel General Monitoring Foam Electrodes**

Recommended for: ICU/CCU, Emergency, O.R., Patient Rooms

This top seller has a fluid-resistant foam backing.

» Flexible foam material hugs body contours, ensuring optimal skin contact for consistent tracings

Item No.	Description	Pkg.
MDSM611503	General Monitoring Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611505	General Monitoring Foam	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611550	General Monitoring Foam	50/pk, 1,000/ct, 6,000/cs



#### **MedGel General Monitoring Cloth Electrodes**

Recommended for: ICU/CCU, Patient Rooms

Provides a gentle fit for patients with sensitive skin.

» Breathable, non-woven cloth conforms to the natural curves of the patient for enhanced comfort

Item No.	Description	Pkg.
MDSM611303	General Monitoring Cloth	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611305	General Monitoring Cloth	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611350	General Monitoring Cloth	50/pk, 1,000/ct, 6,000/cs

#### MedGel General Monitoring Tape Electrodes

Recommended for: ICU/CCU, Emergency, O.R., Patient Rooms

Offers clear visibility of skin under the electrode.

» Transparent, breathable tape lets you evaluate the patient's skin without removing the electrode

Item No.	Description	Pkg.
MDSM611403	General Monitoring Tape	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611405	General Monitoring Tape	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611450	General Monitoring Tape	50/pk, 1,000/ct, 6,000/cs



MDSM611550



MDSM611350



MDSM611450

# **Diagnostic Resting Tabs**

Our resting tab ECG electrodes use a highly conductive gel to ensure clear tracings for an accurate diagnosis.

#### **MedGel Resting Tab ECG Electrodes**

Recommended for: Cardiology, Emergency

Silver/silver chloride adhesive gel delivers excellent adhesion and signal quality.

» Large, solid conductive area allows maximum contact with skin

Item No.	Description	Pkg.
MDSM616101	Resting Tab Electrode, Small	100/pk, 500/bx, 5,000/cs
MDSM616201	Resting Tab Electrode, Adult	100/pk, 500/bx, 5,000/cs



MDSM616101

MDSM616201

### Radiotranslucent

Our carbon snap allows these electrodes to be used during X-rays and MRIs in areas where metal is prohibited.

#### MedGel Radiotranslucent Foam Electrodes

Recommended for: Cath Lab, Radiology

This top seller has a fluid-resistant foam backing.

» Flexible foam material follows natural curves of the patient for maximum skin contact

Item No.	Description	Pkg.
MDSM611803R	Radiotranslucent Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611805R	Radiotranslucent Foam	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611850R	Radiotranslucent Foam	50/pk, 1,000/ct, 6,000/cs



MDSM611850R

#### **MedGel Radiotranslucent Cloth Electrodes**

Recommended for: Cath Lab. Radiology

Offers patients with sensitive skin a comfortable fit.

» Breathable, non-woven cloth gently conforms to the patient's body contours, enhancing comfort

Item No.	Description	Pkg.
MDSM611250R	Radiotranslucent Cloth	50/pk, 1,000/ct, 6,000/cs



MDSM611250R

#### **MedGel Radiotranslucent Tape Electrodes**

Recommended for: Cath Lab, Radiology

Allows clinicians to see the skin without removing the electrode.

» Clear breathable tape lets you evaluate the patient's skin under the electrode

	ı	ı
Item No.	Description	Pkg.
MDSM611650R	Radiotranslucent Tane	50/pk 1000/ct 6 000/cs



MDSM611650R

# Repositionable

These electrodes are designed to adhere firmly even after being repositioned on the patient's skin.

#### **MedGel Repositionable Electrodes**

Recommended for: ICU/CCU, Emergency, O.R., Cath Lab, Radiology, PICU, Patient Rooms

Designed with a conductive hydrogel, our repositionable electrodes are easy to remove and relocate at least once.

- » Provide clear tracings and firm adhesion even after repositioning
- » Carbon snap for standardization throughout your facility
- » Gentle on sensitive skin

Item No.	Description	Pkg.
MDSM611300	Adult, Solid Gel	5/pk, 100/bx, 1,000/cs



MDSM611300

# **Stress/Diaphoretic Foam Electrodes**

Designed specifically for situations requiring a robust adhesive, these electrodes have a comfortable medical-grade foam that resists fluids. The teardrop shape directs sweat around the electrode and stays securely in place during demanding tracing applications.


#### MedGel Stress/Diaphoretic Foam Electrode

Recommended for: Cardiology, Stress Test

Exceptionally powerful adhesive withstands heavy sweating.

» High-performance gel and large conductive area provide clear tracings

Item No.	Description	Pkg.
MDSM611903	ECG Stress Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611905	ECG Stress Foam	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611950	ECG Stress Foam	50/pk, 1,000/ct, 6,000/cs



MDSM611950

#### MedGel Radiotranslucent Stress/Diaphoretic Foam Electrode

Recommended for: Cardiology, Cath Lab, Radiology, Stress Test

Engineered for radiotranslucent situations that require a stronger adhesive.

- » Carbon snap can be used during procedures and in areas where metal is prohibited
- » High-performance gel and extended conductive surface ensure clear tracings

Item No.	Description	Pkg.
MDSM611703R	Radiotranslucent Stress Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611705R	Radiotranslucent Stress Foam	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611750R	Radiotranslucent Stress Foam	50/pk, 1,000/ct, 6,000/cs



MDSM611750R


#### MedGel Wet Gel Stress/Diaphoretic Foam Electrode

Recommended for: Cardiology, Stress Test, EMS

Provides quick reading and stabilization, particularly during heavy patient sweating.

 Highly conductive liquid gel promotes instant, high-quality trace pickup

Item No.	Description	Pkg.
MDSM618303	Wet Gel Stress Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM618305	Wet Gel Stress Foam	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM618350	Wet Gel Stress Foam	50/pk, 1,000/ct, 6,000/cs



MDSM618350

## **Pediatric Electrodes**

Our pediatric electrodes are carefully designed with children in mind. They feature a small footprint, soft materials and greater gel coverage for gentle application and removal.

#### MedGel Pediatric Foam Electrode

Recommended for: PICU

With their smaller size these electrodes are a perfect fit for young patients.

- » High-performance gel and extended conductive surface ensure clear tracings
- » Easy removal with lift tab

Item No.	Description	Pkg.
MDSM611050	Pediatric Foam	50/pk, 1,000/ct, 6,000/cs



MDSM611050

#### MedGel Pediatric Radiotranslucent Electrode

Recommended for: Cath Lab, Radiology, PICU

Our carbon snap allows these electrodes to remain on skin during procedures and in areas where metal is prohibited.

- » Radiotranslucent foam and small size
- » High-performance gel and extended conductive surface provide clear tracings
- » Easy removal with lift tab

Item No.	Description	Pkg.
MDSM6112250R	Pediatric Radiotranslucent Foam	50/pk, 1,000/ct, 6,000/cs



MDSM6112250R

#### MedGel Pre-Wired Neonatal Electrode

Recommended for: NICU

Gentle, breathable tape and ouchless removal with 100% gel coverage make these electrodes perfect for neonatal patients.

- » Three baby-sized electrodes come pre-wired for convenience
- » High-performance gel and extended conductive surface deliver clear tracings every time

Item No.	Description	Pkg.
MDSM61123	Neonatal Clear Tape	3/pk, 300/ct, 3,600/cs
MDSM61143R	Neonatal Radiotranslucent Clear Tape	3/pk, 300/ct, 3,600/cs



MDSM61143R

MedGel Electrodes Recommended Electrodes for Optima This list is not comprehensive and should be used as a recommended guideline only.	al Utilization	ICU/CCU	Emergency	0.R.	Cardiology	Cath Lab	Radiology	Stress Test	PICU	NICU	Patient Rooms	EMS
General Monitoring Foam	+ + + + + + + + + + + + + + + + + + +	•	•	•							•	
General Monitoring Cloth		•									•	
General Monitoring Tape	-	•	•	•							•	
Resting Tab ECG	99		•		•							
Radiotranslucent Foam						•	•					
Radiotranslucent Cloth						•	•					
Radiotranslucent Tape						•	•					
Repositionable	-	•	•	•		•	•		•		•	
Stress/Diaphoretic Foam	<b>(4)</b>				•			•				
Radiotranslucent Stress/Diaphoretic Foam	<b>③</b>				•	•	•	•				
Wet Gel Stress/Diaphoretic Foam	-				•			•				•
Pediatric Foam	<b>*</b>								•			
Radiotranslucent Pediatric Foam	<b>3</b>					•	•		•			
Neonatal Electrodes										•		

For more information and samples, contact your Medline Representative.



Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

Medline Canada 1-800-396-6996 medline.ca | canada@medline.com Medline México

01-800-831-0898 medlinemexico.com | mexico@medline.com

FOLLOW US fin BLOG