


## Intermate System

### Intermate System

The Intermate Systems are portable, disposable devices that use elastomeric technology to infuse drug solutions at a continuous rate. The Large Volume Intermate System (275 mL maximum volume) and the Small Volume Intermate System (105 mL maximum volume) provide a variety of flow rates. Infusion times range from thirty minutes to five hours of continuous infusion. Each unit includes an integral administration set with a Luer lock adapter. Easy to fill; can be used with automated filling systems. Ideally suited for homecare or out-patients requiring intermittent infusion therapy.



### Intermate System – Small Volume

catalogue number	description	maximum volume	flow rate (mL / hr)	nominal delivery time	units / case
2C1734K	Intermate SV 200	105 mL	200 mL / hr	30 minutes	48
2C1732K	Intermate SV 100	105 mL	100 mL / hr	1 hour	48
2C1730K	Intermate SV 50	105 mL	50 mL / hr	2 hours	48

## Intermate System

### Intermate System – Large Volume

catalogue number	description	maximum volume	flow rate (mL / hr)	nominal delivery time	units / case
2C1744K	Intermate LV 250	275 mL	250 mL / hr	1 hour	24
2C1742K	Intermate LV 100	275 mL	100 mL / hr	2.5 hours	24
2C1740K	Intermate LV 50	275 mL	50 mL / hr	5 hours	24

### Intermate System – Extra Large Volume

catalogue number	description	maximum volume	flow rate (mL / hr)	nominal delivery time	units / case
2C1064K	Intermate XLV	550 mL	250 mL / hr	2 hours	12

## Accessories

### Accessories

catalogue number	description	units / case
2C1100	Belt Bag – For use with Small and Large Volume Intermates and Infusors	6