



Flat Skin Barriers


CenterPointLock - Flextend

- ▶ Integrated floating flange
- ▶ Letter matching system

Cut - to - Fit 

		Flange Size	Barrier Opening up to	Product Code	Letter	Quantity per box	Suggested Retail Price
	With Tape Border	1-1/2" (38 mm)	1" (25 mm)	8727	H	5	\$53.84
		1-3/4" (44 mm)	1-1/4" (32 mm)	8722	I		
		2-1/4" (57 mm)	1-3/4" (44 mm)	8723	J		
		2-3/4" (70 mm)	2-1/4" (57 mm)	8724	K		
	Without Tape Border	1-3/4" (44 mm)	1-1/4" (32 mm)	8752	I	5	\$55.03
		2-1/4" (57 mm)	1-3/4" (44 mm)	8753	J		\$76.36
		2-3/4" (70 mm)	2-1/4" (57 mm)	8754	K		
		** No Starting Hole					

Pre - Sized

		Flange Size	Barrier Opening	Product Code	Letter	Quantity per box	Suggested Retail Price
	With Tape Border	1-1/2" (38 mm)	7/8" (22 mm)	8743	H	5	\$57.24
		1-1/2" (38 mm)	1" (25 mm)	8744	H		
		1-3/4" (44 mm)	1-1/8" (29 mm)	8745	I		
		1-3/4" (44 mm)	1-1/4" (32 mm)	8746	I		
		2-1/4" (57 mm)	1-3/8" (35 mm)	8747	J		
		2-1/4" (57 mm)	1-1/2" (38 mm)	8748	J		

Did you know...

- Skin Gel Wipes (7917) or other brands of skin protective wipes should not be used with Flextend Skin Barriers as this may interfere with expected wear times.