

Accessories with Ceramide

Adapt CeraRing Barrier Rings

Flat



Product #	Type	Size	Box of
8805	Standard	2" (48 mm)	10
8815	Slim	2" (48 mm) 2.30 mm thickness	

Convex



Product #	Inner Diameter	Stretch Up To	Box of
89520	13/16" (20 mm)	1" (25 mm)	10
89530	1-3/16" (30 mm)	1-3/8" (35 mm)	
89540	1-9/16" (55 mm)	1-3/4" (45 mm)	

Oval Convex



Product #	Inner Diameter	Stretch Up To	Box of
89601	7/8" x 1-1/2" (22 x 38 mm)	1-1/8" x 1-3/4" (27 x 43 mm)	10
89602	1-3/16" x 1-7/8" (30 x 48 mm)	1-3/8" x 2-1/8" (35 x 53 mm)	
89603	1-1/2" x 2-3/16" (38 x 56 mm)	1-3/4" x 2-3/8" (43 x 61 mm)	

CeraRing Barrier Rings

- Choose round or oval, based on the size and shape of the stoma
- Use convex rings individually or stack them to build a customized level of convexity
- Pair with a flat skin barrier for flexible convexity
- Add to an integrated convex barrier for deeper convexity that may be needed for certain abdomens or stomas