

CURAD® Elastic Net

Compare to Surgilast®*

Rayon/polyester knit. Favored by wound care nurses everywhere. Tubular retainer net holds dressings firmly in place without compression. It is an essential product for fragile, at-risk skin, or for patients with tape allergies.

- » Easy to cut
- » Reduces dressing time
- » Not made with neoprene or natural rubber latex



NONNET04

Item No.	Description	Size	Stretched Width	Relaxed Width	Pkg.
Elastic Net, 25-yd. Stretched Length, Domestic					
NONNET012	Small Fingers, Toes, Wrists, Circumcision	½	4" (10.2 cm)	⅝" (1.6 cm)	1/bx
NONNET01	Large Fingers, Toes, Wrists	1	6" (15.2 cm)	0.6" (1.5 cm)	1/bx
NONNET02	Small Hand, Arm, Leg, Foot	2	7⅔" (19.5 cm)	0.8" (2 cm)	1/bx
NONNET03	Medium Hand, Arm, Leg, Foot, Multiple Fingers	3	9⅓" (23.7 cm)	1" (2.5 cm)	1/bx
NONNET04	Large Hand, Arm, Leg, Foot	4	11⅓" (29.6 cm)	1.2" (3 cm)	1/bx
NONNET05	Average Arm, Small Leg/Knee	5	12⅜" (31.4 cm)	1.5" (3.8 cm)	1/bx
NONNET06	Large Leg, Amputation, Arm Board, Small Head	6	17" (61 cm)	1.9" (4.8 cm)	1/bx
NONNET07	Average Head, Face, Ears, Neck	7	28⅔" (72.8 cm)	2.3" (5.8 cm)	1/bx
NONNET08	Small Chest, Back, Perineum, Axilla, Shoulder, Large Head	8	31⅓" (79.6 cm)	2.8" (7.1 cm)	1/bx
NONNET09	Medium Chest, Back, Perineum, Axilla, Shoulder	9	35" (88.9 cm)	3.2" (8.1 cm)	1/bx
NONNET10	Large Chest, Back, Perineum, Axilla, Shoulder	10	41" (93.1 cm)	3.6" (9.1 cm)	1/bx
NONNET11	X-Large Chest, Back, Perineum, Axilla, Shoulder	11	46" (101.6 cm)	4.3" (10.9 cm)	1/bx
NONNET12	XX-Large Chest, Back, Perineum, Axilla, Shoulder	12	56" (142.2 cm)	4.7" (11.9 cm)	1/bx
NONNET22	Bariatric/Custom Size	22	80" (203.2 cm)	8" (20.3 cm)	1/bx

*This product is not manufactured by Western Medical, Inc., makers of Surgilast®.

TIPS!

If you lose your applicator, simply take a 20cc syringe (without the needle), remove the plunger and place the netting over the syringe. Place the digit inside the syringe and slowly move netting on to the digit.